## **Ophthalmic trial frames**

## **TF-2 Ophthalmic trial frame**

Light-weight trial frame designed to be used with both, adults and children (two different size nose rests included), assures full adjustment against nose and ears position.

PD adjustment range 50-80mm.



## **TF-4 Ophthalmic trial frame**

Solid trial frame designed to be used with adults, assures full adjustment against nose and ears position.

PD adjustment range 48-80mm.



## **TF-2 Ophthalmic trial frame (Pediatric)**

Trial frame with special, adjustable nose pads perfect for examination of children. It's characterised by high durability & excellent quality.

PD adjustment range 48-60mm.



Model	TF-2	TF-4	TF-2 (Pediatric)
Weight (+/-5g)	63g	79g	48g
PD adjustment range (pupil spacing):	50-80mm	48-80mm	48-60mm
Step of the axis scale	5°	5°	5°
Range of lens rotation:	0-180°	0-180°	0-180°
Temple tips angle adjustment range:	+/-10°	+/-10 <sup>0</sup>	+/-15°
Vertex Distance adjustment range:	0-13mm	0-10mm	n/a
Nose pad adjustment range:	0-23mm	0-12mm	0-4mm
Max. quantity of trial lenses	5	4	4

TYPE	DEVICE	PRICE
TF-2	Ophthalmic trial frame	<b>21</b> 8
TF-4	Ophthalmic trial frame	<b>21</b> 8
TF-2 (Pediatric)	Ophthalmic trial frame	<b>23</b> 6